



EU DECLARATION OF CONFORMITY

FOR PERSONAL PROTECTIVE EQUIPMENT

Originator: J.F ROBLES

Revision: 011

Revision date: 20.07.2020

Validity date: 24.02.2025

PRODUCT	SHIELDskin XTREME™ Sterile White Nitrile 400 DI+
DESCRIPTION	Powder-free Extra DI washed Pair-packed Hand-specific 40cm Cleanroom Gloves
CLASSIFICATION	Personal Protective Equipment (PPE) Category III (Complex Design)

SHIELD Scientific codes	Sizes
69 8771	5.5
69 8772	6
69 8773	6.5
69 8774	7
69 8775	7.5
69 8776	8
69 8777	8.5
69 8778	9
69 8779	10

The manufacturer established in the Union:

SHIELD Scientific B.V.

Dr Willem Dreeslaan 1 – 6721 ND BENNEKOM – THE NETHERLANDS

declares under his/her sole responsibility that the PPE (product codes as mentioned above) described hereafter:

SHIELDskin XTREME™ Sterile White Nitrile 400 DI+

is in conformity with the provisions of Regulation (EU) 2016/425 and with the harmonized standards ISO 374-1:2016 (as a Type B glove against reagents: K, P & T), ISO 374-2:2014 (performance level 3, including protection against viruses), EN 16523-1:2015 + A1:2018, ISO 374-4:2013, ISO 374-5:2016 and ISO 21420:2020. This device is identical to the PPE, which is the subject of EU Type Examination (Module B) certificate of conformity *no. F120/81841983* issued by the Notified Body:

SGS FIMKO OY (Notified Body No: 0598)
Takomotie 8, FI-00380 Helsinki, Finland

This device is subject to the procedure set out in Article VIII (Module D) of the Regulation under the surveillance of the Notified Body:

SGS FIMKO OY (Notified Body No: 0598)
Takomotie 8, FI-00380 Helsinki, Finland

Signed for and behalf of SHIELD Scientific B.V



J.F ROBLES
General Manager

Date: 20th July 2020
Place: Bennekom

Validity of this declaration: 20th July 2020 until 24th February 2025